

ELATA<sup>®</sup>

Clinical Foam Mattress E300



Scan for technical resources

# ELATA<sup>®</sup>



## ENGINEERED. TESTED. TRUSTED.

WHERE OTHERS SEE ONE FRAME, WE SEE THE WHOLE PICTURE.

We partnered with **Flinders University** to build a **next generation test system**. Traditional mapping captures one static snapshot every 30 minutes. Our system records **thousands of pressure points** every second. This delivers a **true picture** of how weight shifts and movement affects support.

That depth of data underpins Elata's range of mattress models; E300, E400, E500, E600, and E700. So you can trust the results.

## AUSTRALIAN DESIGNED & ASSEMBLED. LOCALLY TESTED.

**Elata<sup>®</sup> pumps are assembled and rigorously tested here in Australia.** That local focus keeps us agile. We listen, adjust, and roll out updates quickly, strengthening performance and delivering the fast, intelligent care clinicians expect.



**ELATA**

Temperature: 22°C - 71.6°F




Weight: 75 kg

Duration: 60 mins

Proportional Weighted  
95th Percentile Test  
Dummy

## NON-POWERED

- Mattress: E300
- Date: 2025 / 09 / 25
- MAX: 31.86 mmHg
- Average: 3.20 mmHg
- Frame: 600 / 3600



**ELATA<sup>®</sup>**



[support@enablelifecare.com.au](mailto:support@enablelifecare.com.au)

1300 370 370

[enablelifecare.com.au](http://enablelifecare.com.au)